

ORGANOID ISOLATION – large intestine

Crypt Isolation Media DMEM/F12 Advanced
Pen/Strep
Glutamine or Glutamax
HEPES (10mM final)

Post-isolation media (ENYCF) Crypt isolation media +
40 ng/ml Epidermal Growth Factor (Murine Recombinant, Peprotech 315-09),
50nM LDN-193189 (Selleckchem, Catalog No.S2618)
10uM Y27632: ROCK inhibitor (Selleckchem, Catalog No.S1049)
5uM CHIR99021 inhibitor (CAYMAN CHEMICALS - 13122)
10ug/ml FunginTM (InvivoGen, #ant-fn-1)

Crypt Growth Media (WENR) Crypt isolation media +
40 ng/ml Epidermal Growth Factor (Murine Recombinant, Peprotech 315-09),
50nM LDN-193189 (Selleckchem, Catalog No.S2618)
1:20 RSPO1 conditioned media
0.5nM Fc-WNT (ImmunoPrecise, N001 - 100µg/500 µg)

On Ice: Sterile PBS
Sterile Filtered 5mM EDTA/PBS
Matrigel (BD 354234)

- 1) Euthanize mouse, spray (drench) with 70% EtOH and move into the TC hood. Remove entire colon while avoiding the skin and hair, and trim away as much mesentery as possible.
- 2) Open longitudinally and scrape away excess mucous and fecal matter using a glass slide. Place in 20ml of PBS on ice. Shake vigorously to remove any fecal material. Pour off PBS and repeat until the supernatant is clear – approx. 4 washes
- 3) Cut colon into 5-10mm pieces into ice cold EDTA/PBS (10ml) in 50 ml Falcon tube. Place on shaker / roller at 4°C for 60 minutes. In this time prepare falcon tube with filters
- 4) Shake tube vigorously ~10 times and remove supernatant. Replace with ice cold sterile PBS + 200 U DNaseI (20ul/10ml for one fraction) and shake vigorously ~30 times (fraction I). Repeat once (fraction II).
- 5) Filter fractions twice through a 100µm filter into a 50ml conical tube containing 1ml of FBS
- 6) Spin at 1000rpm for 4 minutes.
- 7) Resuspend in ~150µl Crypt Isolation Media and check under the microscope. Crypts appear as thin U-shaped structures. If you have access to Lgr5-EGFP mice, it can be useful to first run the isolation using these animals as crypts are readily identified by GFP staining at the base.
- 8) Add Matrigel to a small volume of crypts – if possible try to keep the volume you add as small as possible (5ul crypts/100ul of Matrigel works well, but you can double that if crypts are sparse). Keep Matrigel on ice at all times.
- 9) Pipette Matrigel droplets into multiple small wells. We commonly use 40ul in each well of a 48 well plate and establish 4-6 wells per isolation – occasionally you will develop fungal or bacterial contamination. Spreading the sample over multiple wells ensures low-level contamination does not affect the entire sample. Let polymerize for 5-10 minutes in the hood, then harden in a 37°C incubator for 15 minutes.
- 10) Add 1ml of EN+Y+Chir+Fungin for 2 days, check the following day for contamination and then after two days swap to WENR media. Replenish WENR media every 2 days.