

**ORGANOID ISOLATION – large intestine**

**Crypt Isolation Media** DMEM/F12 Advanced  
Pen/Strep  
Glutamine or Glutamax  
HEPES (10mM final)

**Crypt Growth Media (ENR50)** DMEM/F12 Advanced containing Pen/Strep, Glutamine, HEPES  
40 ng/ml Epidermal Growth Factor (Murine Recombinant, Peprotech 315-09),  
50 ng/ml Noggin (Murine Recombinant, Peprotech 250-38)  
250ng/ml R-spondin 1 (R&D Systems 3473-RS) or **conditioned media**  
100ng/ml Wnt3a (Millipore GF-160) or **conditioned media**

**For 100ml of media:** 20 ul of EGF  
100 ul of Noggin  
5ml conditioned media  
25ml conditioned media

*NOTE: This recipe uses ~50% the amount of R-spondin and Noggin originally published in Sato et al, 2009. I have maintained wild-type organoids in this media for 6 months without noticeable changes in cell behavior.*

**On Ice:** Sterile PBS  
Sterile Filtered 5mM EDTA/PBS  
Matrigel (BD 354234)

- 1) Euthanize mouse, spray (drench) with 70% EtOH and move into the TC hood. Remove entire colon while avoiding the skin and hair, and trim away as much mesentery as possible.
- 2) Open longitudinally and scrape away excess mucous and fecal matter using a glass slide. Place in 20ml of PBS on ice. Shake vigorously to remove any fecal material. Pour off PBS and repeat until the supernatant is clear – approx. 4 washes
- 3) Cut into ~5-10mm pieces and place in 10ml of ice-cold 5mM EDTA/PBS. Pipette vigorously up and down using 10ml pipette (~15 times). Remove the supernatant.
- 4) Add 10ml of EDTA/PBS, Invert 15 times then agitate gently (rocker, roller, shaker) for 15 minutes at 4C. Remove supernatant and wash once with ice-cold PBS.
- 5) Add 3ml of basal medium containing 500U/ml Collagenase IV. Pipette up and down 5 times using a 5ml pipette and place in a 37C water bath for 30mins
- 6) Crypt Releasing Step:
  - a. Add 10ml PBS and pipette vigorously ~15 times. Collect supernatant containing villi fragments and crypts as Fraction 1.
  - b. Repeat 5a and collect the supernatant as Fraction 2.

*NOTE: In our hands Fraction 2 always contains the best ratio of crypts : villi/debris, but this will depend on multiple factors and should be determined individually for each lab / user. It may be necessary to collect a third or even fourth fraction, depending on initial results.*

- 7) Add 10ml Crypt Isolation Medium containing DNase I (200U), pipette up and down gently to mix and pour (do not pipette) through a 100µm filter. You may have to use 2-3 filters as they may clog depending on amount of debris.
- 8) Filter flow-through on a 70um filter into a 50ml conical tube containing 1ml of FBS.
- 9) Spin at 1000rpm for 4 minutes.

- 10) Resuspend in ~150 $\mu$ l Crypt Isolation Media and check under the microscope.  
Crypts appear as thin U-shaped structures. If you have access to Lgr5-EGFP mice, it can be useful to first run the isolation using these animals as crypts are readily identified by GFP staining at the base.
- 11) Add Matrigel to a small volume of crypts – if possible try to keep the volume you add as small as possible (5ul crypts/100ul of Matrigel works well, but you can double that if crypts are sparse). Keep Matrigel on ice at all times.
- 12) Pipette Matrigel droplets into multiple small wells. We commonly use 40ul in each well of a 48 well plate and establish 4-6 wells per isolation – occasionally you will develop fungal or bacterial contamination. Spreading the sample over multiple wells ensures low-level contamination does not affect the entire sample. Let polymerize for 5-10 minutes in the hood, then harden in a 37<sup>o</sup>C incubator for 15 minutes.
- 13) Add 250 ul of Crypt Growth Media over the gel. Replenish WENR50 media every 2 days.